Malignant pleural mesothelioma (MPM)-specific DNA methylation patterns in patients using liquid biopsies



S. Schmid^{1,2*}, S. Ul Haq^{1*}, L. Zhan¹, M. C. Brown¹, D. Patel¹, F.A. Shepherd¹, N.B. Leighl¹, A.G. Sacher¹, M. de Perrot¹, B.C. J. Cho¹, F. Zaeimi¹, B.H. Lok¹, S.V. Bratman¹, M.S. Tsao¹, M. Früh³, P.A. Bradbury¹, G. Liu¹.

*Contributed equally, 1Princess Margaret Cancer Centre, University of Toronto; 2Inselspital, University Hospital Berne, University of Berne, Switzerland; 3Cantonal Hospital St. Gallen, Switzerland


Background

- The Cancer Genome Atlas (TCGA) show DNA methylation patterns are distinct between tumors and normal tissue
- DNA methylation profiles are highly tissue-specific and can identify tissue-of-origin from cell free DNA (cfDNA).
- Circulating tumor (ct)DNA profiling in MPM has been limited by molecular heterogeneity
- No clear recurrent mesothelioma-specific mutations
- Methylome profiling could identify underlying mesothelioma-specific patterns
- Cell-free methylated DNA immunoprecipitation sequencing (cfMeDIP-seq) can be used on plasma cfDNA
- cfMeDIP-seq on MPM and healthy non-cancer controls (NCC) plasma samples can identify epigenetic changes, overcoming molecular heterogeneity

Methods

- cfMeDIP-seq was applied to pre-treatment plasma samples of 39 MPM patients and 16 NCCs
- cfMeDIP-seq libraries sequenced at a depth of 70 million reads, paired-end, on the NovaSeq 6000.
- For all analyses, chromosomes 1-22 were binned into 300bp windows (total of 9.6e6 genome-wide windows) and reads were tallied per bin
- Differential methylated regions (DMRs) were examined using DESeq2
- Signficant DMRs (log2FC > 1, p-adj < 0.05) were characterized for MPM and NCC samples for CpG and gene body features
- KEGG pathway analysis was also done on signficant DMRs

Results





Baseline characteristics


Volcano plot of significant DMRs in Meso versus NCC plasma samples

Each dot represents a 300bp window. Horizontal line corresponds to p-adj = 0.05. Vertical lines are +/- 1 log2 fold-change.

Yellow: Significant DMRs in Meso. Blue: Significant DMRs in NCC. Grey: Non-significant DMRs







Summary and Conclusion

- Meso samples are enriched for tumor-specific CpG islands and prpmoter regions. KEGG-pathway analysis in MPM samples identified several important pathways. Several promotors significantly hypermethylated in MPM have also been implicated in other cancers
- ctDNA global methylome profiling is a promising tool to identify novel biological pathways and targets in MPM. There is potential for future use of cfMeDIP-Seq in MPM screening, minimal residual disease, and therapeutic monitoring.

Sabine.schmid@insel.ch Abstract #4125